



 GE HealthCare

# Mammography

Imaging Accessories & Supplies



Reliable



Compatible



One-Stop  
Shop

# How to Order



Available 24/7, GE HealthCare Service Shop makes ordering supplies online fast and easy.

- Search real-time inventory of compatible, OEM products
- Take advantage of Auto Reorder for parts you buy regularly
- Batch order multiple items at once with Quick Order
- Get instant order and tracking confirmations
- Pay by PO or credit card or integrate your purchasing system

Visit GE HealthCare Service Shop to register.

1. Visit GE HealthCare Service Shop, navigating to the registration link in the upper left-hand corner of the page.
2. Give us your account details. Pro tip: Provide your customer account number for expedited registration.
3. A welcome message will come to you through the email you provide. Start shopping!

Register now at [ServiceShop.GEHealthcare.com](http://ServiceShop.GEHealthcare.com)

## Other ways to order:


Phone: (866) 281-7545 Options 1, 2, 1, 2

Fax: (877) 279-6990

Email: [gehcaccessorysales@ge.com](mailto:gehcaccessorysales@ge.com)

EDI: Free – contact us for details

Hours: 7AM – 6PM (CST)



# Table of Contents

- Injectors ..... MM 4 ▶
- SensorySuite for Mammography..... MM 6 ▶
- ABUS System Consumables ..... MM 7 ▶
- Biopsy Solutions..... MM 8 ▶
- Mammography Phantom.....MM 11 ▶
- Mammography Furniture.....MM 12 ▶
- Power Solutions .....MM 14 ▶

# Injectors

## Medrad® Stellant FLEX® Injection System for Contrast Enhanced Mammography

The MEDRAD Stellant FLEX helps you meet your financial, clinical, and organizational goals, so that you can continue to deliver quality patient care while maintaining healthy business operations.

### Stellant FLEX Injection System

- Compatible with GE SenoBright HD Contrast Enhanced Spectral Mammography
- User-driven design and easy-to-use interface of MEDRAD Stellant FLEX guides efficient Mammography workflow
- Staff experience with FLEX Contrast Injection System drives satisfaction and standardization

### Stellant FLEX Enhanced Benefits

- Accurate contrast dose administration for better visualization of abnormalities in breast tissue
- Consistent delivery of contrast for enhancement of tumor tissue
- Simple installation to minimize operational disruption
- Small, intuitive and easy to use design to maximize operability

### Certeгра® Workstation

From study set-up and preparation to study administration and results management, the Certeгра Workstation serves as a workflow-centralized technologist interface to help users enhance efficiencies and patient care.



**MEDRAD Stellant FLEX Configuration**

Pedestal-mount configuration includes:

- Dual injector head on pedestal with integral IV pole
- Syringe heat maintainer
- Certegra Workstation with USB drive
- DualFlow software
- Base control unit
- 22.8 m (75 ft) head extension cable
- 7.6 m (25 ft) base to display cable
- Power cord
- Product information package
- Operations manuals
- Installation, customer’s operational training at time of installation, and one year full on-site warranty in Bayer service countries

**More Choices for Operational Efficiency**

Options for Contrast Delivery

- 150 mL syringe option is offered in addition to the 200 mL syringe size for protocols that require higher volumes of contrast
- Syringe options offer you the flexibility to optimize contrast media use, reduce environmental impact, and cut the cost of waste disposal, all of which reduce operational costs
- **Colored beacon technology**  
A clear indication of the presence of fluid in the syringe with an orange colored identifier
- **MEDRAD Fluidots embedded on the syringe**  
Offers additional assurance of fluid in the syringe prior to injection
- **Ultra-clear barrel design**  
Further enhances the visualization of fluid to help prevent user error
- **2D barcodes**  
System uses barcodes to determine: syringe size, expiration date, if the syringe has already been used



| Cat#     | Description  | Warranty |
|----------|--|----------|
| E80141JM | MEDRAD Stellant FLEX Pedestal with Certegra Workstation – For use in Mammography | E        |

# SensorySuite for Mammography

## An interactive experience centered around your patient.



### Sight & Sound

A flat-screen monitor will display soothing video and relaxing ambient sounds of your selected environment.

### Sight

A flat-screen monitor on the wall will display a slide show of photos of your selected ambience.

By focusing on your patients, we've re-invented the experience...and addressed many of the reasons why women avoid mammograms.

Designed to stimulate a woman's senses, SensorySuite can distract her from the perceived discomfort, pain, and anxiety of a mammogram.

SensorySuite lets your patient choose the environmental ambience she prefers for her mammogram:

- Seaside
- Garden
- Waterfall
- Mountain
- Zen Garden
- Natural Rhythms
- Marine
- Mountainside
- Spring Renewal

### Designed to help you stand apart, SensorySuite helps enable you to:

- Improve the mammography experience of every woman
- Increase patient compliance with screening mammography
- Promote the high level of patient care you provide
- Differentiate your healthcare facility from others in your area

Advertise and promote your SensorySuite solution to women and referring physicians in your area with our comprehensive marketing tool kit of promotional materials.

**Take mammography to a whole new place!**

### SensorySuite – Sight & Sound – Interactive Experience

| Cat#     | Description   | Warranty |
|----------|---|----------|
| E63301CJ | Interactive Experience including one small monitor (diag 32 in/81 cm)   | E        |
| E63301CK | Interactive Experience including two small monitors (diag 32 in/81 cm)  | E        |
| E63301CL | Interactive Experience including one large monitor (diag 48 in/122 cm)  | E        |
| E63301CM | Interactive Experience including two large monitors (diag 48 in/122 cm) | E        |

### SensorySuite – Sight – Patient Interactive Device

| Cat#     | Description  | Warranty |
|----------|--|----------|
| E63301CE | Tablet PC for waiting room   | Y        |
| E63301CF | Secured and flexible wall mount system for tablet PC in the waiting room | W        |



E63301CE

# ABUS System Consumables

## ABUS Membranes

The ABUS disposable scan membrane provides comfort to the patient during the 3D image acquisition. The flexible, polyester, mesh fabric stabilizes the breast tissue and enables better image quality during the scan. The membrane is single-use only which eliminates the risk of cross-contamination and the need for sterilization.

### ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo-v™ ABUS Platinum systems

Catalog #E8340AB

Curved membrane dimensions:  
229 mm (L) x 213 mm (W) x 19 mm (D)

## POLYSONIC® Ultrasound Lotion

POLYSONIC® Ultrasound Lotion is an acoustic coupling lotion that is hypo allergenic and has a light pleasing scent.

## Protex™ Ultrasound Disinfectant Wipes

Protex™ one-step cleaner disinfectant kills more than 30 pathogens including HIV, H1N1 and MRSA, while remaining safe for even the most sensitive surfaces. The pre-moistened wipes are easy and convenient to use.



E8340AB



E8365EF

E8365EC

ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo-v™ ABUS (Automated Breast Ultrasound) Platinum systems (non-Sterile)

Compatible for use with either ABUS membrane



E8365EG

E8365EB



E8365EH

Compatible for use with ABUS system

| Cat#    | Description  | Quantity per Box | Warranty |
|---------|--|------------------|----------|
| E8340AB | ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo-v™ ABUS (Automated Breast Ultrasound) Platinum systems (non-Sterile) | 50/box           | W        |
| E8365EC | PolySonic Ultrasound Lotion – (Non-Sterile) 8.5 oz dispenser bottles   | 12 bottles/box   | H        |
| E8365EF | PolySonic Ultrasound Lotion – (Non-Sterile) 1 US gallon with refillable dispenser  | 4 gallons/case   | H        |
| E8365EB | Protex Ultra Disinfectant Wipes – 60 ct Soft-pack w/stay-put adhesive backing – 6.5 in x 6.0 in  | 12 pack/box      | H        |
| E8365EG | Protex Ultra Disinfectant Wipes – 75 ct canister pack - 7 in x 9.5 in  | 8 packs/case     | H        |
| E8365EH | Protex Ultra Disinfectant Wipes – 80 ct canister pack - 7 in x 10 in   | 12/case          | H        |

# Biopsy Solutions

## Ultra DBI 3-in-1 Interventional Mammography Table

The Ultra DBI™ maximizes department space by incorporating upright and decubitus stereotactic breast biopsy, as well as mammography into one bariatric capable platform, increasing throughput and increasing continuity of care.

- Easy transition from chair to table provides an all-encompassing platform for varying patient population
- Accommodates both upright and decubitus positions for stereotactic breast biopsy, as well as mammography and breast ultrasound
- Can also be used as the IB contrast injection chair for CESM (Contrast Enhanced Spectral Mammography) or other interventional procedures
- Optional IV pole accessory for CESM exams
- Bariatric capable patient weight capacity of 500 lbs establishes continuity of care
- Reduce turnover time and increase throughput, even in small exam rooms, with interchangeable head and foot rests eliminating the need to reposition the table to accommodate each patient
- Optional battery for use when a plug-in outlet is not available



E10811SX



E10811SU



E10811SF



E10811SU

| Cat#     | Description  | Warranty |
|----------|--|----------|
| E10811SU | Ultra DBI Table, including Power Cords Kit         | E        |
| E10811SX | Optional battery functionality for Ultra DBI Table | H        |
| E10811SF | IV Pole for Ultra DBI 3-in-1 Mammography Table     | H        |

## Mammography Biopsy Chair

The Hausted® Mammography/Biopsy Chair is designed for use as a procedural chair or treatment table during mammography/biopsy procedures.

### Versatility

- Self storing side rails allow the chair to maneuver into tight spaces
- Small footprint
- Can be used for per-induction, transport, treatment and recovery
- The pneumatic assisted back rest can adjust for supine to 90° for procedural positioning



E6315HD



E6315HD

Seat cushion Width: 22 in/55.88 cm  
 Overall Length (Supine): 75 in/190.5 cm  
 Overall Length (Upright): 42.5 in/107.95 cm  
 Trendelenburg: 10°  
 Maximum Weight Capacity: 325 lbs/147.42 kg

| Cat#    | Description              | Warranty |
|---------|--------------------------|----------|
| E6315HD | Mammography Biopsy Chair | H        |



**Rolling Table Cart**

- May be used in a range of medical applications, including Mammography biopsy procedures
- Wheels allow the cart to be moved from room to room
- Locking casters offer stability when needed
- Height – 35.37 in / 898.5 mm
- Table Top Length – 18.6 in / 473.08 mm
- Table Top Width – 18.3 in / 464.18 mm
- Table Top Depth – 0.20 in / 5.08 mm (this is center/sides have more depth)
- Weight – 114.2 lb / 51.8 kg
- Table Top material – HDPE (High-density Polyethylene)



E6315HD

| Cat#    | Description  | Warranty |
|---------|--|----------|
| E6315SA | Rolling Table Cart for Mammography Biopsy Procedures | C        |

### Sterile Disposable Needle Guides

- Provides a horizontal or lateral approach for core biopsy and localization
- Packaged 2 guides per package, 5 packages per case; sold per case
- Sterile item, disposable – non-returnable and non-refundable



E63101DM

| Cat#                                     | Description   | Needle Size<br>mm | Gauge | Hole Diameter<br>mm | Warranty |
|--|---|-------------------|-------|---------------------|----------|
| <b>Without slot for wire positioning</b> |   |                   |       |                     |          |
| E63101DL                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 21 Gauge, 10/box                                | 0.8               | 21    | 0.9                 | H        |
| E63101DM                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 20 Gauge, 10/box                                | 0.9               | 20    | 1.0                 | H        |
| E63101DN                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 19 Gauge, 10/box                                | 1.1               | 19    | 1.2                 | H        |
| E63101DP                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 18 Gauge, 10/box                                | 1.3               | 18    | 1.4                 | H        |
| E63101DR                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 16 Gauge, 10/box                                | 1.6               | 16    | 1.8                 | H        |
| E63101DS                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 14 Gauge, 10/box                                | 2.0               | 14    | 2.25                | H        |
| <b>With slot for wire positioning</b>    |   |                   |       |                     |          |
| E63101EB                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 20 Gauge with slot for wire positioning, 10/box | 0.9               | 20    | 1.0                 | H        |
| E63101EC                                 | Sterile, disposable STEREOFIX-2 Needle Guide, 19 Gauge with slot for wire positioning, 10/box | 1.1               | 19    | 1.2                 | H        |

### Non-Sterile Reusable Needle Guides

- Non-sterile reusable metal needle guide is to be used for standard mammography biopsy procedures
- Packaged 2 per package, 2 packages per case; sold per case




E63101MC

| Cat#   | Description                                    | Needle Size<br>mm | Gauge | Hole Diameter<br>mm | Warranty |
|--|--|-------------------|-------|---------------------|----------|
| <b>Needle guides with retaining catch</b>    |  |                   |       |                     |          |
| E63101MA                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 0.8               | 21    | 0.9                 | H        |
| E63101MB                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 0.9               | 20    | 1.0                 | H        |
| E63101MC                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 1.07              | 19    | 1.2                 | H        |
| E63101MG                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 2.11              | 14    | 2.25                | H        |
| <b>Needle guides without retaining catch</b> |  |                   |       |                     |          |
| E63101MH                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 2.55              | 13    | 2.55                | H        |
| E63101MJ                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 2.90              | 12    | 2.9                 | H        |
| E63101MK                                     | Non-Sterile Metal Reusable Needle Guide (4/cs) | 3.20              | 11    | 3.2                 | H        |

# Mammography Phantom

## Mammography ACR Phantom

- The phantom simulates the X-ray attenuation of a 4.2 cm slab of compressed human breast composed of 50% adipose tissue and 50% glandular tissue
- An exceptional tool to help identify artifacts due to grids, filters, film-processing (for film-screen systems), as well as those produced by digital mammography systems
- Includes a carrying case



E6322DJ



| Cat#   | Description        | Warranty |
|--|--------------------|----------|
| E6322DJ  | ACR Breast Phantom | H        |
| <b>Compatibility:</b> Analog and digital mammography systems |                    |          |

# Mammography Furniture

## Senographe Pristina™ Accessories Storage Wall Stand

- Holds 8 paddles
- Trays for mag stand and bucky
- Two hooks – for both face shields
- Color coordinated to Senographe Pristina™ system
- Lightweight, aluminum construction for easy wall mounting
- Dimensions – 46.2 (L) x 21 (W) x 17.25 in (D)/1172.9 (L) x 533.4 (W) x 438.9 mm (D)
- Weight – 11.9 lbs

## Senographe Pristina™ Accessories Storage Cabinet

- Slots available for all paddles
- Shelves for mag stands and bucky, or use top of the cabinet for additional storage
- Two hooks – for both face shields
- Color coordinated to Senographe Pristina system
- Sturdy steel base with lightweight aluminum cabinet, shelves and tracks
- Dimensions – 60.2 (L) x 25.3 (W) x 15.4 in (D)/1530.7 (L) x 643.6 (W) x 392.1 mm (D)
- Weight – 57 lb



Both cabinets are specifically designed to store the Senographe Pristina 3D accessories.



E6315TB



E6315TA

| Cat#    | Description  | Warranty |
|---------|--|----------|
| E6315TB | Senographe Pristina Accessories Storage Wall Stand | L        |
| E6315TA | Senographe Pristina Accessories Storage Cabinet    | L        |

**GE Mammography Accessories Cabinet**

Available exclusively from GE HealthCare, this convenient mammography accessory cabinet protects and organizes your mammography components.

- Available for use with SenoClaire™, Senographe Essential™, Senographe DS™, Senographe 2000D and DMR systems
- Holds 14 paddles, including standard or tomo paddles
- Paddles hang vertically to avoid scratching
- Can be wall mounted or free standing
- Cabinet can be used for storing mag stand or phantoms
- Aluminum construction
- Dimensions – 30 in/76.2 cm (L) x 45 in/114.3 cm (H) x 14 in/35.6 cm (D)
- Weight – 49 lb/22.2 kg



Designed for GE tomo paddles



E6315T

| Cat#   | Description                     | Warranty |
|--------|---------------------------------|----------|
| E6315T | Mammography Accessories Cabinet | L        |

**GE Anti-Fatigue Mat**

- Helps to reduce body stress and fatigue resulting from standing for long periods
- One of the most comfortable, buoyant anti-fatigue floor mats on the market
- Durable construction; sealed to prevent moisture absorption and facilitate cleanup



E7058A



E7058AB




E7058AC

| Cat#    | Description   | Warranty |
|---------|---|----------|
| E7058A  | GE Anti-Fatigue Floor Mat (Blue 3' x 5' x 5/8") – 20 lb | H        |
| E7058AB | GE Anti-Fatigue Floor Mat (Gray 3' x 5' x 5/8") – 20 lb | H        |
| E7058AC | GE Anti-Fatigue Floor Mat (Blue 2' x 3' x 5/8") – 9 lb  | H        |

# Uninterruptible Power Supplies

## Uninterruptible Power Supplies (UPS)

- Helps provide reliable, clean, consistent power for Diagnostic Imaging and PACS systems
- Helps to eliminate power anomalies (noise, transients, over and undervoltage), which could damage the system's sensitive computer components
- Enables the DI system to complete scanning after the loss of supply power and allows for the saving of valuable data and orderly system shutdown
- GE HealthCare tests and approves the UPSs listed in this catalog for compatibility with GE HealthCare Diagnostic imaging systems. All solutions are On-Line Double Conversion UPS.
- Configurations are available to support multi-modalities or even complete departments



| Cat#                   | Description  | Warranty |
|------------------------|--|----------|
| <b>Mammography UPS</b> |  |          |
| E4502JJ                | Eaton 6 kVA, 208V, Single Phase UPS (Senographe Pristina only) | 0        |

# Main Disconnect Panels

## Main Disconnect Panels

- Main Disconnect Panels provide emergency shut down, undervoltage protection, overcurrent protection, OSHA lockout tag provisions, and serves as a local disconnect for the GE HealthCare Diagnostic Imaging System.
- Helps reduce installation time by providing a single-point power connection, eliminating the need to mount and wire a number of individual components
- Standardized design and testing helps ensure high product quality and system reliability
- UL and cUL listed for compliance with National Electric Code
- Surface or semi-flush mounting
- Please consult your GE HealthCare Field Engineer to determine if this is the appropriate product for your needs before purchasing



| Cat#                                      | Description   | Warranty |
|---|---|----------|
| <b>Mammography Main Disconnect Panels</b> |   |          |
| E4502PA                                   | Main Disconnect Panel for the Senographe Pristina/Essential/DS (208-240V, 50/60 Hz, UL) | Y        |

## Specifications

- All specifications are subject to change without notice
- The product and brand names herein are the trademarks or registered trademarks of their respective owners

## Pricing

- Please call to obtain current pricing. All pricing and product availability is subject to change without notice
- Price invoiced will be the price in effect at the time of order entry
- Quotations available upon request

## Blanket Purchase Orders

You can save time and money by ordering your consumable products with our Blanket Order program. Blanket orders must be scheduled for a year or more and a signed purchase order agreement must be on file. Contact us for additional details.

## Restocking Fee

- Returns are subject to, but not limited to, a minimum 15% restocking fee
- Certain items, such as sterile and environmentally controlled products within this catalog, are non-returnable and nonrefundable unless defective and have been noted as such

## Product Returns

- Immediately upon receipt, please inspect all deliveries for order accuracy and condition
- For products shipped in error, please obtain return authorization from the Accessories and Supplies Call Center within thirty (30) days of delivery
- Power Quality items are nonreturnable and nonrefundable

## Contact Us for General Information and Warranty Terms

### In U.S.

GE HealthCare  
Accessories and Supplies Customer Support Center  
P.O. Box 414, W-544  
Milwaukee, WI 53202

Phone: (866) 281-7545 – opt 1, opt 2, opt 1, opt 2

Fax: (877) 279-6990

Online: [www.gehealthcare.com/accessories](http://www.gehealthcare.com/accessories)

Email: [gehaccessorysales@ge.com](mailto:gehaccessorysales@ge.com)

### All other regions

Contact your local GE Representative.

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator dedicated to providing integrated solutions to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier.

Learn more at [www.gehealthcare.com](http://www.gehealthcare.com).

GE HealthCare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE HealthCare representative for the most current information and availability. Please visit [www.gehealthcare.com](http://www.gehealthcare.com). Data subject to change.

©2023 General Electric Company.

GE and the GE Monogram are trademarks of General Electric Company.

GE is a trademark of General Electric Company used under trademark license.

All other product names and logos are trademarks or registered trademarks of their respective companies.

JB00642US Rev. 09/23



GE HealthCare